


Product Datasheet

Myosin-10 Antibody (orb1527336)

Description	Myosin-10 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1935 and 1955
Target	Myosin-10
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P35580
NCBI	NP_005955.1
Dilution Range	WB - Not recommended. Use rabbit anti-Myosin-10 antibody A302-043A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human Myosin-10 by western ...