


Product Datasheet

G3BP2 Antibody (orb1527339)

Description	G3BP2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 432 and 482
Target	G3BP2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UN86
NCBI	NP_036429.2
Dilution Range	WB - Not recommended. Use rabbit anti-G3BP2 antibody A302-040A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human G3BP2 by western blot...