

Product Datasheet

ENaC beta antibody (orb152738)

Description

Rabbit polyclonal to ENaC Beta (PerCP). The

membrane ion chann...

Species/Host

Rabbit

Reactivity

Frog, Hamster, Human, Mouse, Rat

Conjugation

PerCP

Tested

ICC, IHC, IP, WB

Applications
Immunogen

Produced against the C-terminal tail (amino acids 617-638) of rat beta ENaC (antibody designation 3755-2)

Target

ENaC

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-404 was sufficient for detection of beta-ENaC in 20 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

87kDa

Uniprot ID
[P37090](#)
NCBI
[NP_036780](#)
Entrez
[24767](#)
Dilution Range

WB (1:1000), IHC (1:100)

Expiration Date

12 months from date of receipt.