


Product Datasheet

TAK1 Antibody (orb1527562)

Description	TAK1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 250 and 300
Target	TAK1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O43318
NCBI	NP_663304.1
Dilution Range	WB - Not recommended. Use rabbit anti-TAK1 antibody A301-916A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TAK1 by western blot ...