



Product Datasheet

ZNF276 Antibody (orb1527626)

Description	ZNF276 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 564 and 614
Target	ZNF276
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8N554
NCBI	NP_001106997.1
Dilution Range	WB - 1:10,000 - 1:25,000; IP - Not recommended
Expiration Date	12 months from date of receipt.



Detection of mouse ZNF276 by western blotting



Detection of human ZNF276 by western blotting