


Product Datasheet

LRWD1 Antibody (orb1527648)

Description	LRWD1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 432 and 482
Target	LRWD1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UFC0
NCBI	NP_690852.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. Use rabbit anti-LRWD1 antibody A301-867A.
Expiration Date	12 months from date of receipt.



Detection of human and mouse LRWD1 by we...