



Product Datasheet

GCN1L1 Antibody (orb1527689)


Description	GCN1L1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 2621 and 2671
Target	GCN1L1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q92616
NCBI	NP_006827.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse and rat GCN1L1 by western blotting.



Detection of human GCN1L1 by western blotting.



Detection of human GCN1L1 by western blotting.