

Product Datasheet

UT-A1 antibody (orb152772)

Description

Rabbit polyclonal to UT-A1 (PerCP). UT-A1, a

the re...

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

PerCP

Tested

ICC, IHC, WB

Applications
Immunogen

Produced against a synthetic peptide mapped to the C-terminal tail (amino acids 911-929) of rat UT-A1 (antibody designation L194)

Target

UT-A1

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-406 was sufficient for detection of UT-A1 in 20 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

97 and 127kDa

Uniprot ID

[Q62668](#)

NCBI

[NP_062220](#)

Entrez

[54302](#)

Dilution Range

WB (1:1000); IHC (1:25)

Expiration Date

12 months from date of receipt.