


Product Datasheet

NF-kappaB1 Antibody (orb1527731)

Description	NF-kappaB1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1 and 50
Target	NF-kappaB1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P19838
NCBI	AAF35232.1
Dilution Range	WB - Not recommended. Use rabbit anti-NF-kappaB1 antibody A301-820A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human NF-kappaB1 by western...