



Product Datasheet

ZCCHC8 Antibody (orb1527755)

Description	ZCCHC8 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 1 and 50
Target	ZCCHC8
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6NZY4
NCBI	NP_060082.2
Dilution Range	WB - Not recommended. Use rabbit anti-ZCCHC8 antibody A301-806A.; IP - 2 - 5 µg/mg lysate; IHC - 1:2,000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.; IHC-IF - 1:1,000 - 1:5,000
Expiration Date	12 months from date of receipt.



Detection of human ZCCHC8 by western blo...



Detection of human ZCCHC8 by immunohisto...