


Product Datasheet

SAS-6 Antibody (orb1527763)

Description	SAS-6 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1 and 50
Target	SAS-6
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6UVJ0
NCBI	NP_919268.1
Dilution Range	WB - Not recommended. Use rabbit anti-SAS-6 antibody A301-802A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human SAS-6 by western blot...