



Product Datasheet

PCAF Antibody (orb1527994)


Description	PCAF Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IP, WB
Immunogen	between 782 and 832
Target	PCAF
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q92831
NCBI	NP_003875.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human and mouse PCAF by wes...



Detection of human PCAF by western blot ...



Localization of PCAF Binding Sites by Ch...