



Product Datasheet

cAbl Antibody (orb1528005)

Description	cAbl Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP
Immunogen	between 850 and 900
Target	cAbl
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P00519
NCBI	NP_005148.2
Dilution Range	WB - Not recommended. Use rabbit anti-cAbl antibody A301-660A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.; ICC-IF - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections. For cytospin preparations of formaldehyde fixed cells permeabilization with Triton-X 100 is recommended.
Expiration Date	12 months from date of receipt.



Detection of human cAbl by western blot ...



Detection of human cAbl by immunohistochem...