


Product Datasheet

SFRS18/SRrp130 Antibody (orb1528086)

Description	SFRS18/SRrp130 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 200 and 250
Target	SFRS18/SRrp130
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8TF01
NCBI	NP_056306.1
Dilution Range	WB - Not recommended. Use rabbit anti-SFRS18/SRrp130 antibody A301-610A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human SFRS18/SRrp130 by wes...