


Product Datasheet

PHC3 Antibody (orb1528160)

Description	PHC3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 933 and 983
Target	PHC3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8NDX5
NCBI	AAI31773.1
Dilution Range	WB - Not recommended. Use rabbit anti-PHC3 antibody A301-569A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human PHC3 by western blot ...