



Product Datasheet

FIP200/RB1CC1 Antibody (orb1528225)

Description	FIP200/RB1CC1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 1544 and 1594
Target	FIP200/RB1CC1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8TDY2
NCBI	NP_055596.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse FIP200/RB1CC1 by west...



Detection of human FIP200/RB1CC1 by WB a...