

Product Datasheet

Amyloid Oligomers (A11) antibody (orb152844)

Description

Rabbit polyclonal to Amyloid Oligomers

proteins can sometimes oligomer...

Species/Host

Rabbit

Reactivity

Human, Mouse, Other, Rat

Conjugation

ATTO-488

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Synthetic molecular mimic of soluble oligomers

Target

Amyloid Oligomers (A11)

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:200 dilution of SPC-506 was sufficient for detection of amyloid oligomers in 10 µg of mouse brain lysates by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

Uniprot ID
P05067
NCBI
NM_000484.2
Entrez
8666
Dilution Range

The suggested dilution for the current lot of SPC-506 in Western Blot and Dot Blot is 1:100

Expiration Date

12 months from date of receipt.