


Product Datasheet

AKAP13/AKAP-Lbc Antibody (orb1528461)

Description	AKAP13/AKAP-Lbc Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 2754 and 2804
Target	AKAP13/AKAP-Lbc
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q12802
NCBI	NP_009131.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human AKAP13/AKAP-Lbc by WB...