


Product Datasheet

AKAP13/AKAP-Lbc Antibody (orb1528463)

Description	AKAP13/AKAP-Lbc Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1850 and 1900
Target	AKAP13/AKAP-Lbc
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q12802
NCBI	NP_009131.2
Dilution Range	WB - Not recommended. Use rabbit anti-AKAP13/AKAP-Lbc antibody A301-404A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human AKAP13/AKAP-Lbc by we...