

Product Datasheet

**Amyloid Oligomers (A11) antibody [Out of stock]
(orb152851)**

Description

Rabbit polyclonal to Amyloid

Amyloid monomeric proteins can
someti...

Species/Host

Rabbit

Conjugation

AP

Tested

ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Synthetic molecular mimic of soluble
oligomers

Target

Amyloid Oligomers (A11)

Preservatives

PBS, 50% glycerol, 0.09% sodium
azide *Storage buffer may change
when conjugated

Storage

Conjugated antibodies should be
stored at 4°C

Note

For research use only

Application notes

A 1:200 dilution of SPC-506 was
sufficient for detection of amyloid
oligomers in 10 µg of mouse brain
lysates by colorimetric immunoblot
analysis using Goat anti-rabbit
IgG:HRP as the secondary antibody.

Clonality

Polyclonal

Uniprot ID
P05067
NCBI
NM_000484.2
Entrez
8666
Dilution Range

The suggested dilution for the
current lot of SPC-506 in Western
Blot and Dot Blot is 1:100

Expiration Date

12 months from date of receipt.