



Product Datasheet

Pericentrin/Kendrin Antibody (orb1528567)


Description	Pericentrin/Kendrin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	between 175 and 225
Target	Pericentrin/Kendrin
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O95613
NCBI	NP_006022.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate; ICC-IF - 1:1,000 - 1:5,000. Formaldehyde or Acetone fixation is recommended. Permeabilization with Triton-X 100 is recommended for formaldehyde-fixed cells.
Expiration Date	12 months from date of receipt.



Detection of human Pericentrin-Kendrin b...



Detection of human Pericentrin-Kendrin b...



ICC: Sample: HeLa cells. Antibody: Rabbi...