

Product Datasheet

Amyloid Fibrils (OC) antibody [Out of stock] (orb152860)

Description	Rabbit polyclonal to Amyloid Fibrils antibody (ATTO-390)
Species/Host	Rabbit
Reactivity	Human
Conjugation	ATTO-390
Tested Applications	DOT, ELISA, ICC, IHC, IP, WB
Immunogen	Fibrils prepared from human amyloid beta 42 peptide
Target	Amyloid Fibrils (OC)
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-507 was sufficient for detection of amyloid fibrils on PVDF membranes using transferred fibrils by colorimetric dot blot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
Dilution Range	WB (1:1000), IHC (1:100)
Expiration Date	12 months from date of receipt.