



Product Datasheet

PAF49 Antibody (orb1528668)

Description	PAF49 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 225 and 275
Target	PAF49
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O15446
NCBI	NP_036231.1
Dilution Range	WB - Not recommended. Use rabbit anti-PAF49 antibody A301-295A.; IP - 2 - 5 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human PAF49 by immunohistoc...



Detection of human PAF49 by western blot...