


Product Datasheet

SRGAP1/ARHGAP13 Antibody (orb1528685)

Description	SRGAP1/ARHGAP13 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 825 and 875
Target	SRGAP1/ARHGAP13
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7Z6B7
NCBI	NP_065813.1
Dilution Range	WB - Not recommended. Use rabbit anti-SRGAP1/ARHGAP13 antibody A301-286A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of
human
SRGAP1/ARHGAP13
by we...