



Product Datasheet

SPTAN1/Alpha II-spectrin Antibody (orb1528757)


Description	SPTAN1/Alpha II-spectrin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 2422 and 2472
Target	SPTAN1/Alpha II-spectrin
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q13813
NCBI	NP_003118.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended; IHC - 1:100 - 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.; IHC-IF - 1:50 - 1:250
Expiration Date	12 months from date of receipt.



Detection of human and mouse SPTAN1/Alph...



Detection of human SPTAN1/Alpha II-Spect...



IHC: Sample: human breast carcinoma. Ant...