



Product Datasheet

SPTAN1/Alpha II-spectrin Antibody (orb1528760)

Description	SPTAN1/Alpha II-spectrin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 1000 and 1050
Target	SPTAN1/Alpha II-spectrin
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q13813
NCBI	NP_003118.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse SPTAN1/Alpha II-Spectrin...



Detection of human SPTAN1/Alpha II-Spectrin...