


Product Datasheet

CHD7 Antibody (orb1528808)

Description	CHD7 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 2025 and 2075
Target	CHD7
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9P2D1
NCBI	NP_060250.2
Dilution Range	WB - Not recommended. Use rabbit anti-CHD7 antibody A301-223A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CHD7 by western blot ...