


Product Datasheet

TRIM3/BERP Antibody (orb1528823)

Description	TRIM3/BERP Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 400 and 450
Target	TRIM3/BERP
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O75382
NCBI	NP_006449.2
Dilution Range	WB - Not recommended. A301-209A has not performed satisfactorily when used for WB of TRIM3/BERP in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of TRIM3/BERP.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TRIM3/BERP by western...