


Product Datasheet

TRIP10/CIP4 Antibody (orb1528869)

Description	TRIP10/CIP4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 25 and 75
Target	TRIP10/CIP4
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q15642
NCBI	AAM46851.1
Dilution Range	WB - Not recommended. Use rabbit anti-TRIP10/CIP4 antibody A301-186A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TRIP10/CIP4 by wester...