


Product Datasheet

PPP2R5D Antibody (orb1529028)

Description	PPP2R5D Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 500 and 550
Target	PPP2R5D
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14738
NCBI	NP_006236.1
Dilution Range	WB - Not recommended. Use rabbit anti-PPP2R5D antibody A301-098A or A301-100A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human PPP2R5D by western blotting