


Product Datasheet

CAF-1 p60 Antibody (orb1529059)

Description	CAF-1 p60 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 375 and 425
Target	CAF-1 p60
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q13112
NCBI	NP_005432.1
Dilution Range	WB - Not recommended. Use rabbit anti-CAF-1 p60 antibody A301-085A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CAF-1 p60 by WB and I...