



Product Datasheet

DDX46 Antibody (orb1529136)

Description	DDX46 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 300 and 350
Target	DDX46
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7L014
NCBI	NP_055644.2
Dilution Range	WB - Not recommended. Use rabbit anti-DDX46 antibody A301-052A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse DDX46 by immunofluorescence.



Detection of human DDX46 by western blotting.