

Product Datasheet

Chicken IgM-Heavy Chain Antibody (orb1529190)

Description	Chicken IgM-Heavy Chain Antibody
Species/Host	Goat
Reactivity	Gallus
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Target	Chicken IgM-Heavy Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:2,000 - 1:20,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:10,000 - 1:100,000
Expiration Date	12 months from date of receipt.