

Product Datasheet

Chicken IgM-Heavy Chain Antibody (orb1529193)

Description	Chicken IgM-Heavy Chain Antibody
Species/Host	Goat
Reactivity	Gallus
Conjugation	AP
Tested Applications	ELISA, ICC, IHC, WB
Target	Chicken IgM-Heavy Chain
Preservatives	50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl ₂ , 0.1 mM ZnCl ₂ containing 0.2% rAlbumin and 0.09% NaN ₃
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB - 1:5,000 - 1:50,000; IHC - 1:100 - 1:1,000; ICC - 1:100 - 1:1,000; ELISA - 1:1,000 - 1:10,000
Expiration Date	12 months from date of receipt.