



Product Datasheet

CBX3 Antibody (orb1529275)

Description	CBX3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 133 and C-term
Target	CBX3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q13185
NCBI	NP_009207.2
Dilution Range	WB - Not recommended; IP - 2 - 5 µg/mg lysate; IHC - 1:2,000 - 10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse CBX3 by imm...



Detection of human CBX3 by western blot ...