



Product Datasheet

CTNND1 Antibody (orb1529289)

Description	CTNND1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 625 and 675
Target	CTNND1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O60716
NCBI	NP_001322.1
Dilution Range	WB - Not recommended. Use rabbit anti-CTNND1 antibody A300-976A or A300-978A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with Tris-EDTA buffer pH 9.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human CTNND1 by immunohisto...



Detection of human CTNND1 by western blo...