



Product Datasheet

ANKRD28/PITK Antibody (orb1529294)

Description	ANKRD28/PITK Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 975 and 1025
Target	ANKRD28/PITK
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O15084
NCBI	NP_056014.2
Dilution Range	WB - 1:2,000 - 1:1,0000; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human ANKRD28/PITK by immun...



Detection of human ANKRD28/PITK by WB an...