



Product Datasheet

FALZ/BPTF Antibody (orb1529295)


Description	FALZ/BPTF Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 1350 and 1400
Target	FALZ/BPTF
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q12830
NCBI	NP_872579.2
Dilution Range	WB - 1:2,000 - 1:1,000; IP - 2 - 10 µg/mg lysate; IHC - 1:200 to 1:1000. Epitope retrieval with Tris-EDTA pH9.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human FALZ/BPTF by western ...



Detection of human FALZ/BPTF by immunohi...



Detection of human FALZ/BPTF by western ...