



## Product Datasheet

### Topo II Beta Antibody (orb1529339)


|                            |                                                                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Topo II Beta Antibody                                                                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                   |
| <b>Tested Applications</b> | IHC, IP, WB                                                                                                                                                    |
| <b>Immunogen</b>           | between 1570 and C-term                                                                                                                                        |
| <b>Target</b>              | Topo II Beta                                                                                                                                                   |
| <b>Preservatives</b>       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide                                                                                           |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                      |
| <b>Concentration</b>       | 200 µg/ml                                                                                                                                                      |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                        |
| <b>Note</b>                | For research use only                                                                                                                                          |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                     |
| <b>Uniprot ID</b>          | <a href="#">Q02880</a>                                                                                                                                         |
| <b>NCBI</b>                | <a href="#">NP_001059.2</a>                                                                                                                                    |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate; IHC - 1:200 to 1:1000. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                |



Detection of human and mouse Topo II Beta...



Detection of human Topo II Beta by immun...



Detection of human Topo II Beta by immun...