



Product Datasheet

Topo II Beta Antibody (orb1529342)


Description	Topo II Beta Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 1475 and 1525
Target	Topo II Beta
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q02880
NCBI	NP_001059.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate; IHC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Topo II Beta by western blotting.



Detection of human Topo II Beta by immunohistochemistry.



Detection of human Topo II Beta by immunoprecipitation and western blotting.