



Product Datasheet

PPP1CA Antibody (orb1529417)

Description	PPP1CA Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 280 and C-term
Target	PPP1CA
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P62136
NCBI	NP_002699.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. Efficiency at IP is very low.; IHC - 1:500 to 1:2,000. Epitope retrieval with Tris-EDTA pH9.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse PPP1CA by w...



Detection of human PPP1CA by immunohisto...