


Product Datasheet

COP1/RFWD2 Antibody (orb1529432)

Description	COP1/RFWD2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 100 and 150
Target	COP1/RFWD2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8NHY2
NCBI	NP_071902.2
Dilution Range	WB - Not recommended. Use rabbit anti-COP1/RFWD2 antibody A300-894A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human COP1/RFWD2 by western...