



Product Datasheet

SUPT5H Antibody (orb1529480)

Description	SUPT5H Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 1040 and C-term
Target	SUPT5H
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O00267
NCBI	NP_003160.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. Use rabbit anti-SUPT5H antibody A300-868A.; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human SUPT5H by western blo...



Detection of human SUPT5H by immunohisto...