



## Product Datasheet

### DHX38 Antibody (orb1529502)

|                            |                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | DHX38 Antibody                                                                                                                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                           |
| <b>Tested Applications</b> | IHC, IP                                                                                                                                                                                                |
| <b>Immunogen</b>           | between 1150 and 1200                                                                                                                                                                                  |
| <b>Target</b>              | DHX38                                                                                                                                                                                                  |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide                                                                                                                                 |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                                                              |
| <b>Concentration</b>       | 1000 µg/ml                                                                                                                                                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                                                  |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q92620</a>                                                                                                                                                                                 |
| <b>NCBI</b>                | <a href="#">NP_054722.2</a>                                                                                                                                                                            |
| <b>Dilution Range</b>      | WB - Not recommended. Use rabbit anti-DHX38 antibody A300-859A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                        |



Detection of human DHX38 by immunohistoc...



Detection of human DHX38 by western blot...