



Product Datasheet

Histone H3 Antibody (orb1529568)


Description	Histone H3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 100 and C-term
Target	Histone H3
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P68431
NCBI	NP_003520.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate; IHC - 1:200 - 1:1,000. Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE tissue sections.; IHC-IF - 1:100 - 1:500
Expiration Date	12 months from date of receipt.



Detection of human Histone H3 by western...



Detection of mouse Histone-H3 by western...



Detection of human Histone H3 by immunoh...