


Product Datasheet

SUPT6H Antibody (orb1529607)

Description	SUPT6H Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1553 and C-term
Target	SUPT6H
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7KZ85
NCBI	NP_003161.2
Dilution Range	WB - Not recommended. Use rabbit anti-SUPT6H Antibodies A300-801A and/or A300-802A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human SUPT6H by western blo...