



Product Datasheet

DHX33 Antibody (orb1529615)

Description	DHX33 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 1 and 50
Target	DHX33
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9H6R0
NCBI	NP_064547.2
Dilution Range	WB - Not recommended. Use rabbit anti-DHX33 antibody A300-800A.; IP - 2 - 10 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human DHX33 by immunohistoc...



Detection of human DHX33 by western blot...