


Product Datasheet

FIGNL1 Antibody (orb1529664)

Description	FIGNL1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 325 and 375
Target	FIGNL1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q6PIW4
NCBI	NP_071399.2
Dilution Range	WB - Not recommended. Use rabbit anti-FIGNL1 antibody A300-772A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human FIGNL1 by western blotting.