



Product Datasheet

GAPDH Antibody (orb1529864)

Description	GAPDH Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 300 and C-term
Target	GAPDH
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P04406
NCBI	NP_002037.2
Dilution Range	WB - 1:2,000 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of mouse GAPDH by western blot...



Detection of human GAPDH by western blot...