



Product Datasheet

SKIV2L2 Antibody (orb1529907)

Description	SKIV2L2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	between 100 and 150
Target	SKIV2L2
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P42285
NCBI	NP_056175.2
Dilution Range	WB - 1:2,000 - 1:10,000; IHC - 1:100 - 1:500. Epitope retrieval with Tris-EDTA pH9.0 is recommended for FFPE tissue sections.; ICC - 1:100 - 1:500. Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE cell sections.
Expiration Date	12 months from date of receipt.



Detection of mouse SKIV2L2 by western bl...



Detection of human SKIV2L2 by western bl...